




Sizing Balloon

FEATURES & BENEFITS



Soft compliance balloon material

- Made of soft compliance TPU plastic membrane
- Suitable for sizing and minimize invasive of defect
- Support low pressure inflation

45° angle to the shaft

- Providing a right-angled position in relation to the atrial septum.


Four Radiopaque marker bands

- Precise measurement : four radiopaque markbands at 10 mm, 5 mm and 2mm intervals
- Used as a distance reference for magnification correction

ORDERING INFORMATION



Code	Balloon Length (mm)	Max Defect Size (mm)	Shaft Size (Fr.)	Usable Length (mm)	Compatible Guidewire (inches)	Max. Volume of Inflation (cc)
LT-SZB-24	40	18	7F	750	0.035	25
LT-SZB-34	40	28	7F	750	0.035	55



For further inquiries, please contact our local sales representatives.


Lifetech Scientific (Shenzhen) Co., Ltd.

Add: Floor 1-5, Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China
 Tel: +86-755-86026250 Fax: +86-755-86026251
 Email: marketing@lifetechmed.com
 Web: www.lifetechmed.com



EU authorized representative:
 LifeTech Scientific (Europe) Coöperatief U.A
 Oliemolenstraat 60, 6416CB Heerlen, the Netherlands

2020-EN-1.0



The AcuMark™ Sizing Balloon is a double-lumen catheter designed to measure cardiovascular defects such as an ASD.

